

## Singleron PythoN™ - Instrument Parameters

Model	SGR-TDAp10
Catalog number	MD1101001
Voltage	220 VAC
Power	250 W
Frequency	50/60 Hz
Operation Mode	Touch Screen
Instrument Size	Length 365 mm, Width 325 mm, Height 360 mm
Instrument Weight	10 kg
Operating Temperature	10-35 °C
Relative Humidity	10-80% Relative Humidity (no condensation)
Air Pressure	700-1060 hpa

## Singleron PythoN™ - Ordering Information

Product Name	Model	Catalog number
Singleron PythoN™ Tissue Dissociator	SGR-TDAp-10	MD1101001

## Related Consumables - Ordering Information

Kit Name	Singleron PythoN™ Tissue Dissociation Kit	
Catalog number	11300601	
Component	sCelLive™ Dissociation Solution	Singleron PythoN™ Dissociation Tubes
Specification	24RXNs	24pcs/box
Package	4 RXNs/bottle; 6 bottles	/
Storage Temp.	-20°C	RT

## Singleron Biotechnologies

**Address:** Gottfried-Hagen-Str. 60, 51105 Cologne, Germany

**Mail:** [info@singleronbio.com](mailto:info@singleronbio.com)

**To learn more, please visit:**

[www.singleronbio.com](http://www.singleronbio.com)

**Phone:**

+49 221 16824777

## Singleron PythoN™ Tissue Dissociation System

Low Input, High Yield, High Viability, Suitable for Biopsies



